

Absaugpumpen - Übersicht

JAXI-05 (MC 600 A):

FEATURES

Single piston cylinder; no lubrication or daily maintenance;
no pollution to the environment

High flow, high vacuum, good performance and
fast vacuum build-up.

The vacuum can be regulated by the vacuum regulator
according to the requirement,

and displayed by the vacuum gauge

Overflow protection device

SPECIFICATIONS

Power source: AC220V \pm 10% /50Hz
Power input : 40 VA
Max. vacuum: \geq 0.08Mpa(600mmHg)
Flow rate: \geq 25L/min
Noise : \leq 55dB(A)
Collection bottle: 1L
GW/NW : 6.0KG/ 4.0KG
Dimensions : 36cm x 18cm x 29cm



JAXI-06 (MC-600C):

Features

High vacuum, high flow,

Oil free vacuum pump AC100V~240V 50/60Hz;

External DC 12V, and internal battery

Overflow protection device Plastic material is handy and beautiful.

The wall-mounting frame makes it widely used in the room,

on vehicle and conveniently hung

on the side of the wheelchair and displayed

Specifications

Power source: AC220V \pm 10%/50Hz, DC12V, or battery
Battery duration
on full charge : \geq 50 min
Power input : 40 VA
Max.Vacuum: \geq 0.08Mpa(600mmHg)
Vacuum range: 0.02Mpa (150mmHg)~Max. vacuum
Flow rate : 25L/min
Noise: \leq 55dB(A)
Collection bottle: 1L
GW/NW : 6.5KG/4.5KG
Dimensions: 36cm \times 18cm \times 29cm

